

Phonak Bolero™ Marvel



Love at first sound

For more information visit www.phonakpro.com

A Sonova brand

PHONAK
life is on

Product description




Bolero M-M
Bolero M-M Trial
M90/M70/M50/M30




Bolero M-PR
Bolero M-PR Trial
M90/M70/M50/M30

Battery size		312			
Multi-function button		•		•	
Telecoil		•		•	
Bluetooth®		•		•	
Nano coating		•		•	
IP rating		IP68 ⁷		IP68 ⁷	
Dimensions	L x W x D	31.3 x 14.1 x 7.3 (1.23 x 0.56 x 0.29")		36.7 x 16.4 x 8.5 (1.44 x 0.65 x 0.33")	
	Weight	2.18 g / 0.077 oz		3.71 g / 0.131 oz	
Charging times ⁸	0% to 100%			3 h	
	0% to 80%			1.5 h	
	30 min. fitting			10 min.	
		SlimTube 4.0	HE11 680	SlimTube 4.0	HE11 680
Max. Power Output (dB SPL)	2 cc coupler	122	128	127	130
Max. gain (dB)	2 cc coupler	58	63	60	66
Frequency range (Hz)		<100 - >6500	<100 - >8000	<100 - >7000	<100 - >6800
Working current (mA)		1.5	1.6		

Fitting range



Mild to moderately-severe hearing loss, all audiometric configurations



Mild to severe hearing loss, all audiometric configurations

■ SlimTube 4.0
■ HE11 680

Housing colours



P1
Sand Beige



P3
Sandalwood



P4
Chestnut



P5
Champagne



P6
Silver Gray



P7
Graphite Gray



P8
Velvet Black



T7
Alpine White



HO
Beige

⁷ IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard.
⁸ Actual charging time depends on the remaining battery life, but it will not exceed 3 hours.

Performance levels

AutoSense OS™ 3.0	Premium M90	Advanced M70	Standard M50	Essential M30
Calm Situation	•	•	•	•
Speech in Noise	•	•	•	•
Comfort in Noise	•	•	•	
Music	•	•		
Speech in Loud Noise ³	•			
Speech in Car	•			
Comfort in Echo	•			
Media speech + mic	•	•	•	•
Media music + mic	•	•	•	•
Additional programs				
Max. additional programs	3	3	3	3
Comfort in Echo	•			
Speech in Loud Noise	•	•		
Speech in 360°	•	•		
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Comfort in Noise	•	•	•	•
Music	•	•	•	•
Acoustic phone	•	•	•	•
Custom program	•	•	•	•
Streaming programs				
RogerDirect + mic	•	•	•	•
PartnerMic + mic	•	•	•	•
Phone call + mic	•	•	•	•

Features	Premium M90	Advanced M70	Standard M50	Essential M30
UltraZoom	•	•	•	•
SNR-Boost	•	•	•	
DuoPhone	•	•	•	
SoundRecover2	•	•	•	•
User Preference Tuning	•	•	•	•
Real Ear Sound	•	•	•	
Fine tuning channels	20	16	12	8
WhistleBlock	•	•	•	•
NoiseBlock	•	•	•	•
WindBlock	•	•		
EchoBlock	•			
SoundRelax	•	•		
QuickSync	•	•	•	•
Acoustically Optimized Vent	•	•	•	•
Tinnitus Balance	•	•	•	•
auto Acclimatization	•	•	•	•
Environmental balance	•	•	•	•
RogerDirect™	•	•	•	•
Roger and Directional setting	•	•	•	

Fitting

	Bolero M-M Bolero M-M Trial M90/M70/M50/M30	Bolero M-PR Bolero M-PR Trial M90/M70/M50/M30
Phonak Target	6.1 or higher	6.1 or higher
Noahlink Wireless	•	•

1 Rodrigues, T. & Liebe, S. (2018). Phonak AutoSense OS™ 3.0. The new and enhanced automatic operating system. Phonak Insight, retrieved from www.phonakpro.com/evidence, accessed October 17th, 2018.

2 Jansen, S. & Woodward, J. (2018). Love at first sound: the new Phonak precalculation. Phonak Insight, retrieved from www.phonakpro.com/evidence, accessed October 17th, 2018.

3 Stewart, E., Rakita, L. & Drexler, J. (2019). StereoZoom Part 1: The benefit of wirelessly connected narrow directionality in Phonak hearing aids for speech intelligibility. Phonak Compendium, retrieved from <http://www.phonakpro.com/evidence>, accessed February 26th, 2019.

4 Winneke, A., Latzel, M. & Appleton-Huber (2018). Less listening- and memory effort in noisy situations with StereoZoom, Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed October 17th, 2018.

5 Legarth, S., Latzel, M. & Rodrigues, T. (2018). Media streaming: The sound quality wearers prefer, Phonak Field Study News, retrieved from www.phonakpro.com/evidence, accessed October 17th, 2018.

6 Thibodeau, L. (2014). Comparison of speech recognition with adaptive digital and FM wireless technology by listeners who use hearing aids. American Journal of Audiology, 23(2), 201-210

Acoustic coupling

Sound outlet

Hook	Damped (HE11 680) transparent
SlimTube 4.0	Available in left and right and five different lengths (00,0,1,2,3)

Ear piece for SlimTube

Dome	Cap dome	
	Open dome	available in three sizes (S, M, L)
	Vented dome	available in three sizes (S, M, L)
	Power dome	available in three sizes (S, M, L)
SlimTip	Acryl	
	Silicone	
	Titanium	

Ear piece for hook

Earmold	
Universal eartip	

Phonak Bolero M Trial™

The new Phonak Bolero M Trial hearing aids can be quickly set to the required performance level in order to suit your recommendation. This provides you the ability to fit your client with the optimal hearing solution already at the first appointment or to provide a replacement when the client's hearing aid is in repair.

Phonak Charger BTE RIC

For charging Bolero Marvel and Audéo Marvel rechargeable hearing aids. The charger features a holding space for large earmolds, slots allowing two hearing aids to be charged simultaneously, a USB-C power inlet and an indicator light.



Phonak Bolero Marvel Bluetooth® versions

Direct calls/audio streaming

Bluetooth® 4.2 wireless technology and most older Bluetooth phones

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.

Wireless communication portfolio

	Bolero M-M Bolero M-M Trial M90/M70/M50/M30	Bolero M-PR Bolero M-PR Trial M90/M70M50/M30
Wireless accessory		
TV Connector	•	•
Phonak PartnerMic™	•	•
Phonak RemoteControl	•	•
Roger™		
Roger™ microphones	•	•
App		
myPhonak	•	•
Phonak Guide	•	•
myCall-to-Text	•	•

For more information on Bolero Marvel please visit www.phonakpro.com/bolero marvel

For information on Roger microphones please visit www.phonakpro.com/roger

For more information on Phonak apps, please visit www.phonakpro.com/smartapps

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonakpro.com